



CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

Non-Radioisotopic Kit for Measuring CaMKK β Activity

CycLex CaMKK beta Kinase Assay Kit

Cat# CY-1185

Intended Use.....	1
Storage.....	1
Introduction.....	2
Principle of the Assay.....	2-3
Materials Provided.....	4
Materials Required but not Provided.....	5
Precautions and Recommendations.....	6
Detailed Protocol.....	7-11
Evaluation of Results.....	12
Assay Characteristics.....	12
Troubleshooting.....	12
Reagent Stability.....	12
Sample Preparation.....	13
Example of Test Result.....	14-16
References.....	17
Related Products.....	17

Intended Use

The CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** is primarily designed to measure the activities of purified Ca²⁺/calmodulin-dependent protein kinase kinase beta (CaMKK β) or recombinant CaMKK β for the rapid and sensitive evaluation of activators or inhibitors. The phospho-threonine specific monoclonal antibody used in this assay kit has been demonstrated to recognize the phospho-threonine 183 in adenosine monophosphate-activated protein kinase (AMPK) α 1/the phospho-threonine 172 in AMPK α 2, which is efficiently phosphorylated by CaMKK β . Additionally, column fractions of cultured primary cells, cell lines, or tissues can be assayed for CaMKK β activity with the CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** if the appropriate dose of CaMKK β specific inhibitor, e.g. STO-609, is used.

Applications of this kit include:

- 1) Screening activators or inhibitors of CaMKK β .
- 2) Evaluating the effects of pharmacological agents on CaMKK β activity *in vitro*.
- 3) Monitoring the purification of CaMKK β activity.

This assay kit is for research use only and not for use in diagnostic or therapeutic procedures.

Storage

- Upon receipt store all components at 4°C.
- Don't expose reagents to excessive light.



Introduction

The AMP-activated protein kinase (AMPK) is a critical regulator of energy homeostasis, and is a potential target for treatment of metabolic diseases as well as cancer (1). Activation of AMPK requires phosphorylation of threonine 172 within the T loop region of the catalytic $\alpha 2$ subunit by the tumor suppressor LKB1 (2, 3) or the Ca^{2+} /calmodulin-dependent protein kinase kinase beta (CaMKK β) (4-6).

CaMKK β can form a complex with and activate AMPK, but CaMKK α cannot. In addition, it was shown that CaMKK β and AMPK associate through their kinase domains, and CaMKK β must be in an active conformation in order to bind AMPK but not to associate with an alternative substrate, Ca^{2+} /calmodulin-dependent protein kinase IV. In contrast to LKB1, the activation of AMPK by CaMKK β does not require an alteration of the ATP:AMP ratio, but rather occurs in response to an increase in intracellular Ca^{2+} . The expression pattern of CaMKK β in cells and tissues is more limited than that of LKB1 and is highest in multiple regions of the brain (7).


Principle of the Assay

The CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** is a semi-quantitative immunoassay for CaMKK β activity. This product can be used to determine the presence of CaMKK β activity in purification column fractions or to follow the kinetics of a purified or partially purified CaMKK β protein as well as screening CaMKK β inhibitor or activator.

The protocol for the quantitative measurement of CaMKK β kinase activity involves incubation of the CaMKK β sample with its substrate, AMPK $\alpha 1$ -GST-fusion protein, in the presence of Mg^{2+} and ATP, followed by transfer this kinase reaction mixture to the well of microtiter plate, which has been pre-coated with a monoclonal antibody specific for phospho-AMPK $\alpha 1$ T183/AMPK $\alpha 2$ T172 for trapping only phosphorylated substrate. The amount of phosphorylated substrate on the well is measured by a horseradish peroxidase conjugated antibody specific for AMPK $\alpha 1$, which then catalyzes the conversion of the chromogenic substrate tetra-methylbenzidine (TMB) from a colorless solution to a blue solution (or yellow after the addition of stopping reagent). The color is quantified by spectrophotometry and reflects the relative amount of CaMKK β activity in the sample. For kinetic analysis, the CaMKK β -containing sample is added to the wells in a similar fashion and at varying times the reaction is stopped by the addition of the chelator, sodium ethylenediaminetetraacetate (EDTA) and the amount of phosphorylated substrate determined as before.


The CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** is designed to accurately determine the presence and relative amount of CaMKK β activity in purification column fractions and to determine non-isotopic kinetic analysis of CaMKK β activity. Careful attention to extraction methods and the assay protocol will provide the investigator with a reliable tool for the evaluation of CaMKK β activity.

Summary of Procedure



Add 100 μ L of reaction mixtures to the non-coated wells.


Incubate for 1 hour at 30°C.



Add of 10 μ L of EDTA Solution to stop kinase reaction.

Transfer 100 μ L of reaction mixtures to the antibody-coated wells.

Incubate for 1 hour at room temp.



Wash the wells.

Add 100 μ L of HRP conjugated anti-AMPK α 1 antibody.

Incubate for 1 hour at room temp.

Wash the wells.

Add 100 μ L of Substrate Reagent.

Add 100 μ L of Stop Solution.

Measure absorbance at 450 nm.



Materials Provided

All samples and standards should be assayed in duplicate. The following components are supplied and are sufficient for the one 96-well microplate kit.

Non-coated Microplate: One microplate supplied ready to use for kinase reaction, with 96 wells (12 strips of 8-wells) in a clear, zip-lock bag.

Antibody-coated Microplate: One microplate supplied ready to use, with 96 wells (12 strips of 8-wells) in a foil, zip-lock bag with a desiccant pack. Wells are coated with anti-phospho-AMPK α 1 T183 / AMPK α 2 T172 monoclonal antibody as a capture antibody.

10X Wash Buffer: One bottle containing 100 mL of 10X buffer containing 2% Tween[®]-20

Kinase Buffer: One bottle containing 20 mL of 1X buffer, used for Reaction Buffer and sample dilution.

20X ATP: One vial of lyophilized ATP Na₂ salt.

20X DTT: Two vials of lyophilized dithiothreitol.

20X BSA: One vial of lyophilized BSA.

100X Calmodulin: One vial of lyophilized calmodulin, used for Reaction Buffer (Ca²⁺/CaM plus).

50X CaCl₂: One vial containing 0.4 mL of 125 mM CaCl₂, used for Reaction Buffer (Ca²⁺/CaM plus).

50X EGTA: One vial containing 0.4 mL of 100 mM EGTA, used for Reaction Buffer (Ca²⁺/CaM minus).

10X AMPK α 1 Substrate: One vial containing 37.5 ng of lyophilized recombinant GST-AMPK α 1 (1-394).

EDTA Solution: One vial containing 2 mL of 0.5 M EDTA, pH 8.0. Ready to use.

HRP conjugated Detection Antibody: One bottle containing 12 mL of HRP (horseradish peroxidase) conjugated anti-AMPK α 1 antibody. Ready to use.

Substrate Reagent: One bottle containing 20 mL of the chromogenic substrate, tetra-methylbenzidine (TMB). Ready to use.

Stop Solution: One bottle containing 20 mL of 1 N H₂SO₄. Ready to use.



Materials Required but not Provided

- **CaMKK beta Positive Control (Full length):** Available from CycLex Co., Ltd., Cat#CY-E1185-2. Unused CaMKK beta Positive Control (Full length) should be stored in aliquots at below -70°C.
- **10X STO-609 (20 µM):** Available from Calbiochem, Cat#570250. 1 mM stock solution (DMSO) diluted 1:50 in Kinase Buffer.
- **Orbital microplate shaker**
- **Pipettors:** 2-20 µL, 20-200 µL and 200-1000 µL precision pipettors with disposable tips
- **Precision repeating pipettor**
- **Wash bottle or multichannel dispenser** for plate washing.
- **Microcentrifuge and tubes** for sample preparation
- **Vortex mixer**
- **Microplate washer:** optional (Manual washing is possible but not preferable)
- **Plate reader:** capable of measuring absorbance in 96-well plates at dual wavelengths of 450 nm/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. The plate can also be read at a single wavelength of 450 nm, which will give a somewhat higher reading.
- **Software package facilitating data generation and analysis:** optional
- **500 or 1000 mL graduated cylinder**
- **Reagent reservoirs**
- **Deionized water of the highest quality**
- **Disposable paper towels**



Precautions and Recommendations

- Allow all the components to come to room temperature before use.
- All microplate strips that are not immediately required should be returned to the zip-lock pouch, which must be carefully resealed to avoid moisture absorption.
- Do not use kit components beyond the indicated kit expiration date.
- Use only the microtiter wells provided with the kit.
- Rinse all detergent residue from glassware.
- Use deionized water of the highest quality.
- Do not mix reagents from different kits.
- The buffers and reagents in this kit may contain preservatives or other chemicals. Care should be taken to avoid direct contact with these reagents.
- Do not mouth pipette or ingest any of the reagents.
- Do not smoke, eat, or drink when performing the assay or in areas where samples or reagents are handled.
- Dispose of tetra-methylbenzidine (TMB) containing solutions in compliance with local regulations.
- Avoid contact with the acidic Stop Solution and Substrate Solution, which contains hydrogen peroxide.
- Wear gloves and eye protection when handling immunodiagnostic materials and samples of rat origin, and these reagents. In case of contact with the Stop Solution and the Substrate Solution, wash skin thoroughly with water and seek medical attention, when necessary.
- **Biological samples may be contaminated with infectious agents. Do not ingest, expose to open wounds or breathe aerosols. Wear protective gloves and dispose of biological samples properly.**
- **CAUTION: Sulfuric Acid is a strong acid. Wear disposable gloves and eye protection when handling Stop Solution.**



Detailed Protocol

The CycLex CaMKK beta Kinase Assay Kit is provided with removable strips of wells so the assay can be carried out on separate occasions using only the number of strips required for the particular determination. Since conditions may vary, running an aliquot of the appropriate positive control (see page 5), separately available from CycLex, should be included in each assay. Disposable pipette tips and reagent troughs should be used for all transfers to avoid cross-contamination of reagents or samples.

Preparation of Working Solution

1. Prepare a working solution of **Wash Buffer** by adding 100 mL of the **10X Wash Buffer** (provided) to 900 mL of deionized (distilled) water (ddH₂O). Mix well. Store at 4°C for two weeks or -20°C for long-term storage.
2. Prepare **20X ATP Solution** by adding **1.6 mL** of **ddH₂O** to the vial of **20X ATP** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **20X ATP Solution** should be **1.25 mM**. Store the solution in small aliquots (e.g. 200 µL) at -20°C.
3. Prepare **20X DTT Solution** by adding **0.5 mL** of **ddH₂O** to the vial of **20X DTT** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **20X DTT Solution** should be **100 mM**. Store the solution in small aliquots (e.g. 100 µL) at -20°C.
4. Prepare **20X BSA Solution** by adding **0.75 mL** of **ddH₂O** to the vial of **20X BSA** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **20X BSA Solution** should be **0.067 mg/mL**. Store the solution in small aliquots (e.g. 100 µL) at -20°C.
5. Prepare **100X Calmodulin Solution** by adding **0.125 mL** of **ddH₂O** to the vial of **100X Calmodulin** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **100X Calmodulin Solution** should be **25 µg/mL**. Store the solution in small aliquots (e.g. 50 µL) at -80°C.
6. Prepare **10X AMPKα1 Substrate Solution** by adding **1.2 mL** of **ddH₂O** to the vial of **10X AMPKα1 Substrate** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **10X AMPKα1 Substrate Solution** should be **31.25 µg/mL**. Store the solution in small aliquots (e.g. 100 µL) at -20°C.
7. Prepare **Positive Control** for kinase reaction.

Dilute the **CaMKK beta Positive Control (Full length)** (CycLex Co., Ltd Cat# CY-E1185-2, see page 5) to the final concentration of ~5 units/µL using **Kinase Buffer** (provided).



8. Prepare **Reaction Buffer** for kinase reaction.

Reaction Buffer (Ca²⁺/CaM plus)

Mix the following reagents.

	96 assays	10 assays	1 assay
Kinase Buffer (provided)	7.2 mL	720 µL	72 µL
10X AMPKα1 Substrate Solution	1.0 mL	100 µL	10 µL
20X ATP Solution	0.5 mL	50 µL	5 µL
20X DTT Solution	0.5 mL	50 µL	5 µL
20X BSA Solution	0.5 mL	50 µL	5 µL
50X CaCl₂ (provided)	0.2 mL	20 µL	2 µL
100X Calmodulin (provided)	0.1 mL	10 µL	1 µL
Total	10 mL	1,000 µL	100 µL

Reaction Buffer (Ca²⁺/CaM minus): for measuring precise CaMKKβ activity (See page 11)

Mix the following reagents.

	96 assays	10 assays	1 assay
Kinase Buffer (provided)	7.2 mL	720 µL	72 µL
10X AMPKα1 Substrate Solution	1.0 mL	100 µL	10 µL
20X ATP Solution	0.5 mL	50 µL	5 µL
20X DTT Solution	0.5 mL	50 µL	5 µL
20X BSA Solution	0.5 mL	50 µL	5 µL
50X EGTA (provided)	0.2 mL	20 µL	2 µL
ddH₂O	0.1 mL	10 µL	1 µL
Total	10 mL	1,000 µL	100 µL

Reaction Buffer (Ca²⁺/CaM/ATP minus): for measuring precise CaMKKβ activity (See page 11)

Mix the following reagents.

	96 assays	10 assays	1 assay
Kinase Buffer (provided)	7.2 mL	720 µL	72 µL
10X AMPKα1 Substrate Solution	1.0 mL	100 µL	10 µL
20X DTT Solution	0.5 mL	50 µL	5 µL
20X BSA Solution	0.5 mL	50 µL	5 µL
50X EGTA (provided)	0.2 mL	20 µL	2 µL
ddH₂O	0.6 mL	60 µL	6 µL
Total	10 mL	1,000 µL	100 µL

80-90 µL of Reaction Buffers per assay well will be needed. Mix well.
Discard any unused Reaction Buffers after use.



Standard Assay

1. Remove the appropriate number of wells of **Non-coated Microplate** and **Antibody-coated Microplate** from the pouch and place them into the well holder. Return any unused wells to the foil pouch, refold, seal with tape and store at 4°C.
2. Prepare all samples (diluted with Kinase Buffer as needed). All samples should be assayed in duplicate.
3. Add **10 µL** of **sample** and **Positive Control** (See “7. Prepare Positive Control for phosphorylation.” on page 7) to each well of **Non-coated Microplate** on ice. The **Positive Control** should be included in duplicate wells in each assay.
4. Start the kinase reaction by addition of **90 µL Reaction buffer (Ca²⁺/CaM plus)** to each well of the **Non-coated Microplate**, cover with plate sealer, and incubate at 30°C for 60 minutes shaking at ca. 300 rpm on an orbital microplate shaker.
5. Stop the kinase reaction by addition of **10 µL** of **EDTA Solution** to each well.
6. Transfer **100 µL** of the reaction mixture in each well of the **Non-coated Microplate** to each well of the **Antibody-coated Microplate**, cover with plate sealer, and incubate at room temperature (ca.25°C) for 60 minutes shaking at ca. 300 rpm on an orbital microplate shaker.
7. Wash wells five times with Wash Buffer making sure each well is filled completely. Remove residual Wash Buffer by gentle tapping or aspiration.
8. Pipette **100 µL** of **HRP-conjugated Detection Antibody** to each well, cover with the lid, and incubate at room temperature (ca.25°C) for 60 minutes shaking at ca. 300 rpm on an orbital microplate shaker.
9. Wash wells five times as same as in step 7.
10. Add **100 µL** of **Substrate Reagent** to each well and incubate at room temperature (ca.25°C) for 5–20 minutes shaking at ca. 300 rpm on an orbital microplate shaker.
11. Add **100 µL** of **Stop Solution** to each well in the same order as the previously added Substrate Reagent.
12. Measure absorbance in each well using a spectrophotometric plate reader at dual wavelengths of 450/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. Read the plate at 450 nm if only a single wavelength can be used. Wells must be read within 30 minutes of adding the Stop Solution.



CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

Note-1: Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.

Note-2: Reliable signals are obtained when either O.D. values do not exceed 0.25 units for the blank (no enzyme control), or 2.5 units for Positive Control.

Note-3: If the microplate reader is not capable of reading absorbance greater than the absorbance of Positive Control, perform a second reading at 405 nm. A new O.D. values, measured at 405 nm, is used to determine Positive Control of off-scale samples. The readings at 405 nm should not replace the on-scale readings at 450 nm.

Recommendations

Special considerations for screening inhibitors

In order to estimate the inhibitory effect on individual CaMKK β activity in the test chemicals correctly, it is necessary to conduct the control experiment of "Vehicle control" at least once for every experiment and "Inhibitor control" at least once for the first experiment, in addition to "Test sample", as indicated in the following table. When test chemicals cause an inhibitory effect on CaMKK β activity, the level of A450 is weakened as compared with "Vehicle control".

For inhibitor screening

Assay reagents	Test sample for inhibitor	Vehicle control	Inhibitor control
Reaction buffer (Ca ²⁺ /CaM plus)	80 μ L	80 μ L	80 μ L
10X Inhibitor or equivalent	10 μ L	-	-
Vehicle for inhibitor	-	10 μ L	-
10X STO-609 (20 μ M) *	-	-	10 μ L
Positive Control ** or CaMKK β sample	10 μ L	10 μ L	10 μ L

* See "Materials Required but not Provided" on page 5.

** 5 units/ μ L of CaMKK beta Positive Control (Full length): See "Materials Required but not Provided" on page 5 and "7. Prepare Positive Control for kinase reaction." on page 7.

1. Following the above table, add the Reagents to each well of the microplate. Finally, initiate reaction by adding 10 μ L of "Positive Control" or "CaMKK β sample" to each well and mixing thoroughly at room temperature. Cover with plate sealer or lid, and incubate at 30°C for 60 minutes shaking at ca. 300 rpm on an orbital microplate shaker.

2. Follow the steps 5-12 of "Standard Assay" on page 9.

Note: Although we suggest to conduct experiments as outlined in the table above, the optimal



CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

experimental conditions will vary depending on the parameters being investigated, and must be determined by the individual user. Especially, an appropriate amount of the enzyme must be optimized by titration of the enzyme and setting the amount, which shows OD value does not exceed plateau range in dose response curve.

Special considerations for measuring precise CaMKKβ activity

In order to measure the activity of CaMKKβ correctly, it is necessary to conduct the control experiment of “Inhibitor control” at least once for every experiment, and “Ca²⁺/CaM minus control” and “ATP minus control” at least once for the first experiment, in addition to “No enzyme control” as indicated in the following table. Although the level of A450 increases in “Test sample” when CaMKKβ enzyme activity is in the sample, the high level of A450 is not observed in “Inhibitor control”, “ATP minus control” and “No enzyme control”.

Assay reagents	Test sample	Inhibitor control	Ca ²⁺ /CaM minus control	ATP minus control	Positive control	No enzyme control
Reaction buffer (Ca ²⁺ /CaM plus)	80 μL	80 μL	-	-	80 μL	80 μL
Reaction buffer (Ca ²⁺ /CaM minus)	-	-	80 μL	-	-	-
Reaction buffer (Ca ²⁺ /CaM/ATP minus)	-	-	-	80 μL	-	-
10X STO-609 (20 μM) *	-	10 μL	-	-	-	-
Vehicle for 10X STO-609	10 μL	-	10 μL	10 μL-	10 μL-	10 μL
CaMKKβ sample	10 μL	10 μL	10 μL	10 μL	-	-
Positive Control **	-	-	-	-	10 μL	-
Buffer for CaMKKβ sample	-	-	-	-	-	10 μL

* See “Materials Required but not Provided” on page 5.

** 5 units/μL of CaMKK beta Positive Control (Full length): See “Materials Required but not Provided” on page 5 and “7. Prepare Positive Control for kinase reaction.” on page 7.

- Following the above table, add the Reagents to each well of **Non-coated Microplate**. Finally, initiate the reaction by adding 10 μL of “**CaMKKβ sample**” or “**Positive Control**” or “**Buffer for CaMKKβ sample**” to each well and mixing thoroughly at room temperature. Cover with plate lid, and incubate at 30°C for 60 minutes shaking at ca. 300 rpm on an orbital microplate shaker.
- Follow the steps 5-12 of “Standard Assay” on page 9.

Note: Although we suggest to conduct experiments as outlined in the table above, the optimal experimental conditions will vary depending on the parameters being investigated, and must be determined by the individual user. Especially, an appropriate amount of the enzyme must be optimized by titration of the enzyme and setting the amount, which shows OD value does not exceed plateau range in dose response curve.



Evaluation of Results

Average the absorbance values for the Positive Control and all experimental sample duplicate values (when applicable). When Positive Control (See “7. Prepare Positive Control for kinase reaction.” on page 7) is included in 50 units/assay as an internal control for the kinase reaction, the absorbance value should be greater than 1.0 with a background less than 0.25 when using Reaction buffer.

Assay Characteristics

The CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** has been shown to detect the activity of purified CaMKK β and column fractions containing CaMKK β . The assay may be used to follow the purification of CaMKK β or may be used to detect the presence of CaMKK β in cell lysates.

Troubleshooting

1. The Positive Control should be run in duplicate, when a standard assay is being performed, using the protocol described in the “**Detailed Protocol**”. Incubation times or temperatures significantly different from those specified may give erroneous results.
2. The reaction curve is nearly a straight line if the kinetics of the assay is of the first order. Variations in the protocol can lead to non-linearity of the curve, as can assay kinetics of other than first order. For a non-linear curve, point to point or quadratic curve fit methods should be used.
3. Poor duplicates, accompanied by elevated values for wells containing no sample, indicate insufficient washing. If all instructions in the “**Detailed Protocol**” were followed accurately, such results indicate a need for washer maintenance.
4. Overall low signal may indicate that desiccation of the plate has occurred between the final wash and addition of Substrate Reagent. Do not allow the plate to dry out. Add Substrate Reagent immediately after wash.

Reagent Stability

All of the reagents included in the CycLex Research Product **CycLex CaMKK beta Kinase Assay Kit** have been tested for stability. Reagents should not be used beyond the stated expiration date. Upon receipt kit reagents should be stored at 4°C. Antibody-coated Microplate should be stored in the original foil bag sealed by the zip lock and containing a desiccant pack.

For research use only, not for use in diagnostic or therapeutic procedures



Sample Preparation

Numerous extraction and purification methods can be used to isolate CaMKK β . The following protocols have been shown to work with a number of different tissues and enzyme sources and are provided as examples of suitable methods. Crude samples can frequently be used without dilution while more concentrated or highly purified CaMKK β should be diluted. It is strongly advised that the user always perform an initial experiment to determine the proper dilution to be used in subsequent experiments. This need not be any more than a single time point assay using serial dilutions of the crude extract, cell lysate or sample fraction taken prior to a purification step. One eight well strip of the plate should be sufficient for this initial experiment. All sample preparation should be performed at 4°C and recovered fractions should be kept at 4°C to prevent loss of enzymatic activity.

CAUTION: It should be noted that this assay kit detects not only CaMKK β activity but also other protein kinases in crude extract and column sample. The presence of CaMKK β protein in the samples should be traced by other methods, e.g. western blotting.

Immunoprecipitation Protocol Followed by Measuring CaMKK β Activity

Preparation of Solution and Reagent

A. Preparation of Cell Lysis Buffer

20 mM Tris HCl, pH 7.5, 250 mM NaCl, 10 % glycerol, 0.5 % Nonidet[®] P-40, 1 mM EDTA, 1 mM EGTA, 0.2 mM PMSF, 1 μ g/mL pepstatin, 0.5 μ g/mL leupeptin, 5 mM NaF, 2 mM Na₃VO₄, 2 mM β -glycerophosphate, 1mM DTT

B. Preparation of Protein A Agarose Beads

Add 5 mL of 1X PBS to 1.5 g of protein A agarose beads. Shake 2 hours at 4°C; spin down. Wash beads twice with PBS. Resuspend beads in 1 volume of PBS. (Can be stored for 2 weeks at 4°C)

Preparing Cell Lysates

1. Plate adherent cells in 10cm dish plate at $\sim 3.0 \times 10^6$ cells/plate and incubate the plate at 37°C for 12-16 hours in CO₂ incubator.
2. Remove media, and wash cells with ice-cold PBS and aspirate.
3. Add 1 mL of ice-cold Cell Lysis Buffer to the plate and incubate on ice for 5 minutes.
4. Scrape off and transfer the lysate to a microcentrifuge tube.
5. Rotate the tube for 60 minutes at 4°C and microcentrifuge at 15, 000 rpm for 10 minutes at 4°C.
6. Transfer the supernatant to a new tube. The supernatant is the cell lysate. If necessary, the lysate can be stored at -70°C.

Immunoprecipitation

1. Take 100 μ L cell lysate and add anti-CaMKK β antibody (Abcam, Cat#ab168818; 1-2 μ g) incubate with gentle rocking for 2 hours or overnight at 4°C.



CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

2. Add protein A agarose beads (20 μ L of 50% bead slurry). Incubate with gentle rocking for 1-3 hours at 4°C.
3. Microcentrifuge for 30 seconds at 4°C. Wash the protein A agarose beads 3 times with 500 μ L of Cell Lysis Buffer and successively once with Kinase Buffer. Keep on ice during washes.
4. Resuspend the protein A agarose beads with 20-40 μ L of Kinase Buffer and use 10 μ L as an enzyme sample* to measure CaMKK β activity according to the procedure in the **step 3-12** of “**Standard Assay**” on page 9.

* *Please take care to transfer the protein A agarose beads to the well of the **Antibody-coated Microplate** as little as possible.*



Example of Test Results

Fig.1 Dose dependency of recombinant CaMKK β




Fig.2 Km for ATP of recombinant CaMKK β






Fig.3 Effect of NaCl on the kinase activity of recombinant CaMKK β





Fig.4 Effect of CaMKK β inhibitor, STO-609, on the kinase activity of recombinant CaMKK β






CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

Fig.5 Effect of EGTA and Ca²⁺/ Calmodulin on native CaMKKβ immunoprecipitated from C2C12 cell





CaMKK beta Kinase Assay Kit

User's Manual

For Research Use Only, Not for use in diagnostic procedures

References

1. Hawley SA, Pan DA, Mustard KJ, Ross L, Bain J, Edelman AM, Frenguelli BG, Hardie DG. Cell Metab. 2005; 2(1): 9-19
2. Shaw RJ, Kosmatka M, Bardeesy N, Hurley RL, Witters LA, DePinho RA, Cantley LC. PNAS USA. 2004; 101(10): 3329-3335.
3. Hawley S, Boudeau J, Reid J, Mustard K, Udd L, Makela T, Alessi D, Hardie DG. J Biol. 2003; 2(4): 28.
4. Hurley RL, Anderson KA, Franzone JM, Kemp BE, Means AR, Witters LA. JBC. 2005; 280(32): 29060-29066.
5. Hawley SA, Pan DA, Mustard KJ, Ross L, Bain J, Edelman AM, Frenguelli BG, Hardie DG. Cell Metab. 2005; 2(1): 9-19.
6. Woods A, Dickerson K, Heath R, Hong S-P, Momcilovic M, Johnstone SR, Carlson M, Carling D. Cell Metab. 2005; 2(1): 21-33.
7. Anderson KA, Means RL, Huang Q-H, Kemp BE, Goldstein EG, Selbert MA, Edelman AM, Freneau RT, Means AR. JBC 1998; 273(48): 31880-31889.

Related Products

- * CycLex CaMKK beta Kinase Assay Kit: Cat#CY-1185
- * CaMKK beta Positive Control (Full length): Cat# CY-E1185-2

PRODUCED BY

CycLex Co., Ltd.
1063-103 Terasawaoka
Ina, Nagano 396-0002
Japan
Fax: +81-265-76-7618
e-mail: info@cyclex.co.jp
URL: <http://www.cyclex.co.jp>

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from CycLex Co., Ltd.. To inquire about licensing for such commercial use, please contact us via email.