



## CycLex's Protocol P-002

### Isolation of crude nuclei from rat liver

1. Wash 50 gram of liver once with cold PBS.
2. Cut the liver into small pieces (c.a. 3 mm cubic) with scissors on ice.
3. Add 100 ml of ice-cold homogenize buffer (20 mM potassium phosphate, pH 7.2, 0.25 M sucrose, 4 mM CaCl<sub>2</sub> 4 mM beta-mercaptoethanol).
4. Homogenize in a Dounce homogenizer in ice (usually 15-20 up and down strokes).
5. Layer 20 ml of the homogenate on 10 ml of sucrose cushion (0.88 M sucrose, 3 mM MgCl<sub>2</sub>, 4 mM CaCl<sub>2</sub>).
6. c.f.g. at 3,000 rpm (1,610 x g) for 10 min at 4 °C.
7. Take a pellet of crude nuclei.
8. Resuspend the pellet in 10 vols. of homogenize buffer (20 mM potassium phosphate, pH 7.2, 0.25 M sucrose, 4 mM CaCl<sub>2</sub> 4 mM beta-mercaptoethanol).
9. Layer the suspension over an equal volume of 0.88 M sucrose solution.
10. c.f.g. at 3,000 rpm (1,610 x g) for 10 min at 4 °C to obtain a pellet of clear nuclei.