



## Extraction of Crude Nuclear Extract for HDAC assay

### Buffers

#### \*Lysis Buffer

- 10 mM Tris HCl pH7.5
- 10 mM NaCl
- 15 mM MgCl<sub>2</sub>
- 250 mM Sucrose
- 0.5 % NP-40
- 0.1 mM EGTA

#### \*Sucrose cushion

- 30 % Sucrose
- 10 mM Tris HCl pH7.5
- 10 mM NaCl
- 3 mM MgCl<sub>2</sub>

#### \*Extraction buffer

- 50 mM HEPES KOH, pH 7.5, 420 mM NaCl, 0.5 mM EDTA Na<sub>2</sub>, 0.1 mM EGTA, 10 % glycerol.

### Procedure

#### Isolation of Nuclei

1. Suspend  $1 \times 10^7$  cells (100 mm dish sub-confluent) into 1ml of lysis buffer.
2. Vortex for 10 second.
3. Keep on ice for 15 min.
4. Spin the cells through 4 ml of sucrose cushion at 1,300 x g for 10 min at 4 C.
5. Discard the supernatant.
6. Wash the nuclei pellet once with cold 10 mM Tris HCl pH7.5, 10 mM NaCl.

#### Extraction of Nuclei

1. Suspend the isolated nuclei in 50-100 ul of extraction buffer.
2. Sonicate for 30 seconds.
3. Stand on ice for 30 min.
4. c.f.g. 15,000 rpm for 10 min.
5. Take supernatant (the crude nuclear extract).
6. Determine protein conc. by Bradford method or equivalent.
7. Store the crude nuclear extract at  $-70$  C until use.

**Note: Do not use any kind of protease inhibitor!!**