



LMW-PTP/ACP1 Positive Control

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures

LMW-PTP/ACP1 Positive Control

(Human, full length, recombinant enzyme expressed in *E. coli*.)

Cat# CY-E1358

Lot No.
For 100 Assays
(20 units/ μ L x 500 μ L)

Product Description: Human full length LMW-PTP/ACP1, containing a *N*-terminal GST tag, expressed in *E. coli*. Purified by GSH agarose chromatography. The LMW-PTP/ACP1 Positive control is designed to use for CycLex Protein Phosphatase LMW-PTP/ACP1 Fluorometric Assay Kit [Cat# CY-1358]. The LMW-PTP/ACP1 Positive Control should be added to the well at 100 units/well. Unused LMW-PTP/ACP1 Positive control should be stored at -70°C. AVOID FREEZE/THAW CYCLES!

Formulation: The LMW-PTP/ACP1 Positive Control is supplied frozen in a buffer containing 50mM Tris-HCl (pH 7.0), 2 mM DTT and 50% glycerol.

Source: Human full length LMW-PTP/ACP1, containing *N*-terminal GST tag, expressed in *E. coli*. The GenBank Accession No. is BC106011.

Molecular Weight: recombinant LMW-PTP/ACP1 demonstrates a single 44 kDa band by SDS-PAGE analysis.

Purity: LMW-PTP/ACP1 is greater than 80 % pure as determined by SDS-PAGE analysis.

Substrates: LMW-PTP/ACP1 dephosphorylates a number of substrates, including phosphorylated PDGF-R, EphA2, EphB2 and synthesized phosphatase substrates.

Inhibitors: Sodium Orthovanadate are known as effective PTP inhibitor.

Specific Activity: 20 units/ μ L. This unit value is determined at the point of production and may vary with time and various conditions. Specific Activity also varies among production lots.

Unit Definitions: One unit is defined as the amount of phosphatase required to release 1 pmol of phosphate from synthesized phosphatase substrate, OMFP (3-o-methylfluorescein phosphate), per minute at 30°C.

Storage: Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, AVOID REPEATED HANDLING AND MULTIPLE FREEZE/THAW CYCLES.

Stability: Unopened vial at -70 °C, for 1 year after delivery.



LMW-PTP/ACP1 Positive Control

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures

Related Products:

- * PTP1B Fluorometric Assay Kit : Cat# CY-1350
- * TC-PTP Fluorometric Assay Kit : Cat# CY-1351
- * Protein Phosphatase Cdc25A Fluorometric Assay Kit: Cat# CY-1352
- * Protein Phosphatase Cdc25B Fluorometric Assay Kit: Cat# CY-1353
- * Protein Phosphatase Cdc25C Fluorometric Assay Kit: Cat# CY-1354
- * Protein Phosphatase Cdc25 Combo Fluorometric Assay Kit: Cat# CY-1355
- * Protein Phosphatase Cdi1/KAP Fluorometric Assay Kit: Cat# CY-1356
- * Protein Phosphatase LMW-PTP/ACP1 Fluorometric Assay Kit: Cat# CY-1358
- * Recombinant PTP1B : Cat# CY-E1350
- * Recombinant TC-PTP : Cat# CY-E1351
- * Recombinant Cdc25A (Catalytic domain): Cat# CY-E1352
- * Recombinant Cdc25B (Catalytic domain): Cat# CY-E1353
- * Recombinant Cdc25C (Catalytic domain): Cat# CY-E1354
- * Recombinant Cdi1/KAP: Cat# CY-E1356
- * Recombinant LMW-PTP/ACP1 : Cat# CY-E1358

General References:

1. Dissing J, Johnsen AH, Sensabaugh GF. Human red cell acid phosphatase (ACP1). The amino acid sequence of the two isozymes Bf and Bs encoded by the ACP1*B allele. J Biol Chem. 1991 Nov 5;266(31):20619-25.
2. Wo YY, McCormack AL, Shabanowitz J, Hunt DF, Davis JP, Mitchell GL, Van Etten RL. Sequencing, cloning, and expression of human red cell-type acid phosphatase, a cytoplasmic phosphotyrosyl protein phosphatase. J Biol Chem. 1992 May 25;267(15):10856-65.

PRODUCED BY

CycLex Co., Ltd.
1063-103 Terasawaoka
Ina, Nagano 396-0002
Japan
Fax: +81-265-76-7618
e-mail: info@cyclex.co.jp
URL: <http://www.cyclex.co.jp>

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from CycLex Co., Ltd. To inquire about licensing for such commercial use, please contact us via email.