



Baculovirus envelope gp64 Mouse Monoclonal Antibody

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures

Baculovirus envelope gp64 Mouse Monoclonal Antibody (Clone AF-2C8)

Cat# CY-M1026

50 µg (1 mg/ml x 50 µL)

Clone Name	Applications	Species Cross-reactivity	Molecular Wt.	Source Isotype
AF-2C8	IC, F, E	N/A	64 kDa	Mouse IgG2a

Background: Recombinant baculoviruses derived from the Autographa californica nuclear polyhedrosis virus (AcNPV) are widely used to express heterologous genes in insect cells, but the use of the baculovirus expression vector system is hampered by slow and tedious procedures for the selection and propagation of baculovirus and for titer determination. Titration of baculoviral stocks is important because when doing experiments, it is important to have consistency between samples, and to achieve the right level of transient expression and when producing viral stocks it is important to know the titer of infectious particles for successful virus production. This antibody can be used for determining titers of baculovirus stocks by immunocytochemical or immunofluorescence method.

Specificity/Sensitivity: Baculovirus envelope gp64 antibody detects the gp64 envelope protein of baculovirus Autographa californica nuclear polyhedrosis virus (AcNPV).

Source/Purification: Monoclonal antibody is produced by immunizing mice with wild type AcNPV particles. IgG is purified by protein A-sepharose chromatography.

Recommended Antibody Dilutions: Immunofluorescence assay for detection of AcNPV-based recombinant baculovirus: 0.5-1 µg/mL

Storage: Supplied in 20 mM phosphate buffer (pH 7.5), 300 mM NaCl, 50 % glycerol. Store at -20°C.

Applications Key: WB:Western Blotting IP:Immunoprecipitation IHC:Immunohistochemistry IC:Immunocytochemistry

F:Flow cytometry E:ELISA FP:Fluorescence Polarization assay

Species Cross-Reactivity Key: H:human M:mouse R:rat Hm:hamster Mk:monkey Mi:mink C:chicken X:*Xenopus* Z:zebra fish All:all species expected Species enclosed in parentheses are predicted to react based on 100% sequence homology. N/A: Not Applicable



Baculovirus envelope gp64 Mouse Monoclonal Antibody

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures

Selected References:

1. Hohmann, A. W. and P. Faulkner. 1983. Monoclonal antibodies to baculovirus structural proteins: determination of specificities by Western blot analysis. *Virology* 125(2): 432-44.
2. Volkman, L. E. and P. A. Goldsmith. 1988. Resistance of the 64K protein of budded *Autographa californica* nuclear polyhedrosis virus to functional inactivation by proteolysis. *Virology* 166(1): 285-9.
3. Blissard, G. W. and G. F. Rohrmann 1989. Location, sequence, transcriptional mapping, and temporal expression of the gp64 envelope glycoprotein gene of the *Orgyia pseudotsugata* multicapsid nuclear polyhedrosis virus. *Virology* 170(2): 537-55.
4. Plonsky, I., M. S. Cho, et al. (1999). An analysis of the role of the target membrane on the Gp64-induced fusion pore. *Virology* 253(1): 65-76.

Fig.1 Immunofluorescence detection of baculovirus infected cells





Fig.2 Sandwich ELISA for measurement of recombinant AcMNPV-based recombinant baculovirus





Baculovirus envelope gp64 Mouse Monoclonal Antibody
Product Data Sheet
For Research Use Only, Not for use in diagnostic procedures

Related Products

* Anti-Adenovirus Hexon: Cat# CY-M1027

PRODUCED BY

CycLex Co., Ltd.
1063-103 Terasawaoka
Ina, Nagano 396-0002
Japan
Fax: +81-265-76-7618
e-mail: info@cyclex.co.jp
URL: <http://www.cyclex.co.jp>

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from CycLex Co., Ltd.. To inquire about licensing for such commercial use, please contact us via email.